

# Anti-H7(A/Netherland/03(H7N7))



Product Name	Anti- H7(A/Netherland/03(H7N7))
Catalog No	IT-003-008M8
Source	Mouse monoclonal antibody, clone 1D2
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA (A/Netherlands/219/03(H7N7)) protein (Genebank No. <a href="#">AAR02640</a> ).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H7 (A/Netherland/03(H7N7)) and H7 (A/England/96 (H7N7)) proteins as tested so far. Minimal cross-reactivity to other subtype.
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-003-008M8**

