

Anti-H7(A/Netherland/03(H7N7))



Product Name	Anti- H7(A/Netherland/03(H7N7))
Catalog No	IT-003-008M9
Source	Mouse monoclonal antibody, clone 22B8
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA (A/Netherlands/219/03(H7N7)) protein (Genebank No. AAR02640).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H7 (A/Netherland/03(H7N7)), H7 (A/England/96 (H7N7)), H15 (A/Australia/83/00(H15N8)); and weakly reacts with H9(A/Hong Kong/1073/99(H9N2)) proteins as tested so far. Minimal cross-reactivity to other subtype.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-008M9

