




Anti-HA(H9N2)(Avian)

Product Name	Anti-HA(H9N2)(Avian)
Catalog No	IT-003-009
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hemagglutinin (amino acid 1-345)(H9N2)(A/chicken/HongKong/NT366/03 (H9N2)) (Genebank accession #AAU11163)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian H9 (H9N2). Cross-reactivity to other HAs not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing H9

Note: