

Anti- HA(H9N2)



Product Name	Anti- HA(H9N2)
Catalog No	IT-003-009M1
Source	Mouse monoclonal antibody
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA1 protein (H9N2)(A/Guinea fowl/Hong Kong/WF10/99 (H9N2)) protein (a.a.20-339) (Genebank accession AAO46082).
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts with most HA(H9N2) strains. No cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.