

Anti-H9(A/Hong Kong/33982/2009)



Product Name	Anti- H9(A/Hong Kong/33982/2009)
Catalog No	IT-003-009M3
Source	Mouse monoclonal antibody, clone 4C2
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA (A/Hong Kong/33982/2009) protein(a.a. 17-339)(Genebank #: EPI241740)
Applications	Neut., ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts H9(A/Hong Kong/33982/2009), and weakly reacts with H9(A/Hong Kong/1073/99(H9N2)) proteins as tested so far. Minimal cross-reactivity to other subtype.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

