



HA₁(H₁N₁)(Influenza A)

Product Name	HA ₁ (H ₁ N ₁) (A/Brisbane/59/2007)
Catalog No	IT-003-0012p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA ₁ (H ₁ N ₁) (A/Brisbane/59/2007) protein (a.a.18-343) (Genebank No. ACA28846)
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(H₁N₁) protein

Note: