




HA₁ (H1N1)

Product Name	HA ₁ (H1N1)
Catalog No	IT-003-001p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS only or in PBS with 0.1% BSA and 25% glycerol.
Protein	6xHis tagged hemagglutinin (amino acid 18-345) (H1N1) (A/New Caledonia/20/99(H1N1)) (Genebank No. AAP34324)
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified HA₁ (H1N1) protein.