




HA₁ (H7N7)

Product Name	HA ₁ (H7N7)(A/chicken/Netherlands/1/03(H7N7))	
Catalog No	IT-003-008p	
Source	Viral protein purified from 293 cell culture	
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.	
Protein	6xHis tagged HA ₁ (H7N7)(A/chicken/Netherlands/1/03(H7N7)) protein (a.a.26-344) (Genebank accession # AAR02639)	
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.	
Purity	≥95%	
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.	
Use Limitation	For research use only, not for use in diagnostic procedures.	
<div><div><div><div>250 kD—</div><div>150 kD—</div><div>100 kD—</div><div>75 kD—</div><div>50 kD—</div><div>37 kD—</div><div>25 kD—</div><div>20 kD—</div><div>15 kD—</div><div>10 kD—</div></div><div><div>E</div></div><div>E: purified HA₁ (H7N7) protein</div></div></div>		