



## HA( $\Delta$ TM)(A/Hong Kong/2108/03(H9N2))

---

Product Name	HA( $\Delta$ TM)(A/Hong Kong/2108/03 (H9N2))
Catalog No	IT-003-0093 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM)( A/Hong Kong/2108/03 (H9N2))protein (a.a.18-522) (Genebank No. ABB58945)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM)(A/Hong Kong/2108/03(H9N2))