



# HA<sub>1</sub> (H9N2)

Product Name	HA <sub>1</sub> (H9N2)(A/guinea fowl/Hong Kong/WF10/99(H9N2))	
Catalog No	IT-003-009p	
Source	Viral protein purified from 293 cell culture	
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.	
Protein	6xHis tagged HA <sub>1</sub> (H9N2)(A/guinea fowl/Hong Kong/WF10/99 (H9N2)) protein (a.a.20-339) (Genebank accession # AAO46082)	
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.	
Purity	≥ 95%	
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.	
Use Limitation	For research use only, not for use in diagnostic procedures.	
	<div> <div> 250 kD— 150 kD— 100 kD— 75 kD— 50 kD— 37 kD— 25 kD— 20 kD— 15 kD— 10 kD— </div> <div> <p>E</p> </div> </div>	E: purified HA <sub>1</sub> (H9N2) protein