



Anti-Human RBP-4

Product Name	Anti- Human RBP-4
Catalog No	IT-010-0011
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Peptide from human retinol binding protein 4(hRBP4) (a.a. 203-212) (CDGRSERNLL).
Applications	WB(1:1000), ELISA(1:2000-1:10000), IP, etc.
Specificity	Reacts with peptide of human RBP-4 (a.a. 203-212) (CDGRSERNLL).
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.

ELISA Profile of IT-010-0011

