

Anti-GP(Marburg/Musoke)



Product Name	Anti-GP(Marburg/Musoke)
Catalog No	IT-014-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing ectodomain of GP (Marburg/Musoke) protein (a.a.24-649) (Genebank No. DQ217792).
Applications	ELISA(1:1000-1:100,000), WB(1:200-1:2000), IP, etc.
Specificity	Reacts with most GP(Marburg) glycoproteins. Not cross-react with GP(Ebola).
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

