



Anti-NP(Ebola/Zaire)

Product Name	Anti-NP(Ebola/Zaire)
Catalog No	IT-014-006
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing the nucleoprotein of (Ebola/Zaire) (a.a.1-739) (Genebank No. AF086833).
Applications	ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with NP (Ebola/Zaire). Cross-react with NP from other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

