

# Anti-G1 (Rift Valley Fever Virus)



Product Name	Anti- G1 (Rift Valley Fever Virus)
Catalog No	IT-017-001M1
Source	Mouse monoclonal antibody, clone 1D8
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Glycoprotein (RVFV) (a.a. 154-584 of the M protein. This region is also called Phlebovirus_G1) (Genebank Accession No: AAA47449).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with glycoprotein (RVFV). Cross-reactivity not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-017-001M1

