



Anti-NP(Rift Valley Fever Virus)

Product Name	Anti-NP (Rift Valley Fever Virus)
Catalog No	IT-017-002M1
Source	Mouse monoclonal antibody
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant NP (Rift Valley Fever Virus) (Strain ZH-548 M12) (a.a. 1-245) (Genebank No. P21700)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with NP (Rift Valley Fever Virus). Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

