

# Monkey IgG



---

Product Name	Monkey IgG ( <i>Macaca mulatta</i> )
Catalog No	IT-000-004
Source	Rhesus monkey (SIV, HBV, HIV etc free)
Isotype	IgG
Formulation	Each vial contains 1 mg purified IgG in 1 ml PBS
Applications	Negative control for ELISA, Neut., IP, IF etc.
Storage	Store at 4°C; stable for 3 months from the date of shipment. For long term storage, aliquot and freeze at -80°C. <b>Avoid repeated freeze and thaw.</b> Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.