

Anti-*Gaussia* Luciferase



Product Name	Anti- <i>Gaussia</i> Luciferase
Catalog No	IT-000-014M2
Source	Mouse monoclonal antibody
Isotype	IgM
Formulation	Each vial contains 100 µg purified IgM in 100 µl PBS with 40% glycerol
Immunogen	Recombinant <i>Gaussia</i> Luciferase (a.a. 1-168) (Genebank No. AAG54095)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with <i>Gaussia</i> Luciferase. No Cross-reactivity to other luciferase protein.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-000-014M2

