

Pre-M (WNV)



Product Name	Pre-M (WNV)
Catalog No	IT-006-052Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in 20mM phosphate buffer pH 7.5. Reconstitute the protein with 100 µl of Millipore water.
Protein	The E.Coli derived 20kda recombinant protein contains the West-Nile N-Terminal Pre-M Virus immunodominant regions. The protein is fused with 6xHis tag.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.