

# Anti-His tag



---

Product Name	Anti-His tag
Catalog No	IT-000-His
Source	Mouse monoclonal antibody
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant protein containing 6xHis tag
Applications	ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with 6xHis tag in tagged proteins
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.