

Anti- 6xHis



Product Name	Anti- 6xHis tag
Catalog No	IT-000-His-Biotin
Source	Mouse monoclonal antibody
Isotype	IgG1
Formulation	Each vial contains 50 µg biotinylated IgG in 200 µl PBS with 40% glycerol
Immunogen	Recombinant protein containing 6xHis tag
Applications	ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with 6xHis tag in tagged proteins
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

