




Anti-NA (H₇N₇)

Product Name	Anti-NA(H ₇ N ₇)
Catalog No	IT-003-028
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	NA from influenza A virus (H ₇ N ₇) (A/chicken/Netherlands/1/03 (H7N7)) (Genebank accession #AAR11365)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with influenza A NA (H ₇ N ₇) protein. Cross-reactivity to other NA proteins not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing NA (H₇N₇)

Note: