

Anti-N7(H7N7)



| | |
|----------------|--|
| Product Name | Anti-N2 (H7N7) |
| Catalog No | IT-003-028M2 |
| Source | Mouse monoclonal antibody, clone 43E4 |
| Isotype | IgG1 |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol |
| Immunogen | Recombinant NA protein (A/chicken/Netherlands/1/03(H7N7)) (Genebank accession No. AAR11365) |
| Applications | ELISA, IP, IF etc. |
| Epitope | Conformational |
| Specificity | Reacts with most contemporary N7 (H7N7). Cross-reacts with N2 (H3N2). No cross-reactivity to other subtypes. |
| Storage | Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

