




Anti-NP (Influenza B)

Product Name	Anti-NP(Influenza B)
Catalog No	IT-003-031
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	NP form influenza B virus (B/Malaysia/2506/2004)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with NP(Influenza B) protein. Cross-reactivity to other clades not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing NP (Influenza B)

Note: