




Anti-HA (Influenza B)

Product Name	Anti-HA(Influenza B)
Catalog No	IT-003-032
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HA from influenza B virus (B/Malaysia/2506/2004)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with Influenza B HA protein. Cross-reactivity not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing HA (Influenza B)

Note: