

Anti- HA (B/Malaysia/2506/2004)



Product Name	Anti- HA(B/Malaysia/2506/2004)
Catalog No	IT-003-032M1
Source	Mouse monoclonal antibody
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA (B/Malaysia/2506/2004)(GeneBank No. ABU99194) protein (a.a. 12-547)
Applications	ELISA, Neut. etc.
Epitope	Conformational
Specificity	Reacts with HA (B/Malaysia/2506/2004)
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.