

Anti-HA(Influenza B(Seasonal))



Product Name	Anti- HA(Influenza B(Seasonal))
Catalog No	IT-003-032M2
Source	Mouse monoclonal antibody, clone 4F12
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HAΔTM (B/Brisbane/60/2008)(Genebank No. ACN29380)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with HA(B/Brisbane/60/2008) and HA(B/Malaysia/2506/2004). Cross-reactivity with other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-032M2

