

Anti-HA(B/Wisconsin/01/2010)



Product Name	Anti-HA(B/Wisconsin/01/2010)
Catalog No	IT-003-032M7
Source	Mouse monoclonal antibody, clone 20B8
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA (B/Wisconsin/01/2010) (a.a. 15-548) (Genebank #: EPI271545).
Applications	Neut., ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with HA(B/Wisconsin/01/2010) only as tested so far. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-032M7

