

# Anti-gH/gL/gp42 Complex(EBV/B95-8)



Product Name	Anti- gH/gL/gp42 Complex(EBV/B95-8)
Catalog No	IT-005-030M11
Source	Mouse monoclonal antibody, clone 8F3
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Co-immunization with DNA vaccine expressing the ectodomains of gH protein (a.a.1-529) (Genebank #: P03231), gL protein (1-137) (Genebank #: CAD53428), and gp42 protein (a.a. 1-223) (Genebank #: YP_401672)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with gH/gL/gp42 complex(EBV/B95-8). Not reactive to gH/gL complex (EBV/B95-8).
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-005-030M11**

